STOTZ FEMORAL SHAFT REVISION CENTRALIZER

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
MEDICAL DEVICE FUTURE SOLUTIONS & DEVELOPMENTS

STOTZ CEMENT EXTRACTION SET

SET REF.: 44-910-00

Femoral Shaft Revision Centralizer
Set, consisting of:

44-912-25	Guide wire with drill tip Ø 2.5mm x length 565mm (PACK 5 PCS.)
44-910-08	Centralizer Ø 8.0mm, cannulated
44-910-10	Centralizer Ø 10.0mm, cannulated
44-910-12	Centralizer Ø 12.0mm, cannulated
44-910-14	Centralizer Ø 14.0mm, cannulated
60-014-06	Medullary reamer, rigid with trinkle connection Ø 8.0mm
003-0200-010TR-AE	Medullary reamer, flexible with trinkle connection Ø 8.0mm
55-198-05	Tray, steel
55-831-50	Silicone slide, fixation inlays, 2 PCS.
C-111-100	1/1 Size container bottom and lid perforated, silver outside 580x280x100 mm inside 550x265x84 mm weight 3156 g

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The evolution of revision:

For bone retention and reduction of trauma for patients

The flexible shaft of the centring aid enables individual alignment of the centring head and therefore an optimum contact position of the drill tip on the cement seal.




DRILLING Safe and minimally invasive

A guide wire with a self-centring drill tip is threaded through the centring instrument and fixed at an optimum placement position at the cement seal.

The guide wire fixed by guide shaft and centring tip can therefore be centrally inserted into the cement seal – this prevents unwanted insertion of the drill tip into the cortical bone.

After removal of the centring aid, the cannulated, flexible drill shaft included in the set with a front-cutting tip is guided into the seal by the prepared guide wire.

The resulting drill channel now allows the cement filling to be removed by conventional intermedullary reamers.

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
INNOVATION

Stotz Femoral Shaft Revision Centralizer – optimises hip-joint revision

This newly designed instrument set manufactured from biocompatible implant steel simplifies targeted removal of the seal of the prosthesis tip without opening the cortical bone of the femur in hip-joint revisions. This prevents injury to the femoral bone and reduces surgical and rehabilitation times and the risk of infection.

Tivolilaan 205
6824 BV Arnhem
The Netherlands
T: +31 -(0)24-3618804
F: +31 -(0)26-3882550
E: info@hb-medical.nl
I: www.hb-medical.nl