



ARTHROSCOPY







Our **DR SHAVER BLADES & BURS**
are set apart by:

- Very good cutting properties
- No metal abrasion
- Long life time and durability
- Fitting hubs
- More convenience for operating personnel
- More precision in surgical manipulations
- Savings on disposable items
- Shorter intervention times
- Reusability

DR SHAVER BUR CERATITAN®



DR BARREL BUR®



DR ROUND BUR®

Internal tube coated with CERATITAN® (an anti-bacterial ceramic coating)

- CERATITAN® - coating for highest
 - cut quality
 - cutting performance
 - contrast by visualization
- with 2 aspiration windows (prevents "snowstorm effect")

DR SHAVER BLADE pat.

- DR WAVE CUTTER® with 1 double edge
- DR PIT BULL® with 3 double edges
- High cut quality and cutting performance
- Suitable for soft tissue and bone resection (PIT BULL®)
- Milling recommended up to 5000 rpm (maximum 8000 rpm)
- Oscillations up to 5000 rpm
- Either clockwise or anti-clockwise rotation
- Available in straight or 15° curved geometry
- Various lengths: 130mm, 155mm and 174mm
- Available single use and reusable
- Supplied sterile



DR WAVE CUTTER®



DR PIT BULL® CERATITAN®



DR PIT BULL®



DR SYNOVIA PIT BULL®

DR WAVE CUTTER (1 double edge)

Diam. Ø (mm)	Allowed rpm	Allowed oscillation rpm	APPLICATIONS			
			knee	shoulder	ankle	hip
4	5000	5000	x	x		
5	5000	5000	x	x		x

DR PIT BULL (3 double edges)

Diam. Ø (mm)	Allowed rpm	Allowed oscillation rpm	APPLICATIONS			
			knee	shoulder	ankle	hip
3	5000	5000	x	x	x	
4	5000	5000	x	x		
5	5000	5000	x	x		x
7	5000	5000				x

DR BARREL BUR

Diam. Ø (mm)	Allowed rpm	Allowed oscillation rpm	APPLICATIONS			
			knee	shoulder	ankle	hip
4	5000	0	x	x	x	x
5.5	5000	0		x		x

DR ROUND BUR

Diam. Ø (mm)	Allowed rpm	Allowed oscillation rpm	APPLICATIONS			
			knee	shoulder	ankle	hip
4	5000	0	x		x	x
5.5	5000	0				x



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES

DROR

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	027.101.230	DROR blue
DR WAVE CUTTER	Ø 5.0	130	027.201.230	DROR blue
DR PIT BULL	Ø 3.0	130	027.301.240	DROR blue
DR PIT BULL	Ø 4.0	130	027.101.240	DROR blue
DR PIT BULL	Ø 5.0	130	027.201.240	DROR blue
DR PIT BULL	Ø 7.0	130	027.701.240	DROR blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.236	DROR blue
DR BENDING PIT BULL	Ø 3.0	130	027.301.340	DROR blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.241	DROR blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.242	DROR blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.250	DROR blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.251	DROR blue
DR BARREL BUR	Ø 4.0	130	137.100.220	DROR blue
DR BARREL BUR	Ø 5.5	130	137.200.220	DROR blue
DR ROUND BUR	Ø 4.0	130	127.100.220	DROR blue
DR ROUND BUR	Ø 5.5	130	127.200.220	DROR blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	027.101.030	DROR
DR WAVE CUTTER	Ø 5.0	130	027.201.030	DROR
DR PIT BULL	Ø 3.0	130	027.301.040	DROR
DR PIT BULL	Ø 4.0	130	027.101.040	DROR
DR PIT BULL	Ø 5.0	130	027.201.040	DROR
DR PIT BULL	Ø 7.0	130	027.701.040	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.036	DROR
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	027.201.043	DROR
DR PIT BULL	Ø 4.0	155	027.101.037	DROR
DR PIT BULL	Ø 5.0	155	027.201.037	DROR
DR BENDING PIT BULL	Ø 3.0	130	027.301.140	DROR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.041	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.042	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.050	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.051	DROR
DR BARREL BUR	Ø 4.0	130	137.100.020	DROR
DR BARREL BUR	Ø 5.5	130	137.200.020	DROR
DR ROUND BUR	Ø 4.0	130	127.100.020	DROR
DR ROUND BUR	Ø 5.5	130	127.200.020	DROR

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	027.101.528	DROR
DR WAVE CUTTER	Ø 5.0	130	027.201.528	DROR
DR PIT BULL	Ø 4.0	130	027.101.538	DROR
DR PIT BULL	Ø 5.0	130	027.201.538	DROR
DR PIT BULL	Ø 7.0	130	027.701.538	DROR
DR PIT BULL	Ø 4.0	155	027.101.535	DROR
DR PIT BULL	Ø 5.0	155	027.201.535	DROR
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.534	DROR
DR BARREL BUR	Ø 4.0	130	137.100.519	DROR
DR BARREL BUR	Ø 5.5	130	137.200.519	DROR
DR ROUND BUR	Ø 4.0	130	127.100.519	DROR
DR ROUND BUR	Ø 5.5	130	127.200.519	DROR

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES

FMS

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	020.101.230	FMS blue
DR WAVE CUTTER	Ø 5.0	130	020.201.230	FMS blue
DR PIT BULL	Ø 3.0	130	020.301.240	FMS blue
DR PIT BULL	Ø 4.0	130	020.101.240	FMS blue
DR PIT BULL	Ø 5.0	130	020.201.240	FMS blue
DR PIT BULL	Ø 7.0	130	020.701.240	FMS blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	020.101.236	FMS blue
DR BENDING PIT BULL	Ø 3.0	130	020.301.340	FMS blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	020.101.241	FMS blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	020.101.242	FMS blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	020.101.250	FMS blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	020.101.251	FMS blue
DR BARREL BUR	Ø 4.0	130	130.100.220	FMS blue
DR BARREL BUR	Ø 5.5	130	130.200.220	FMS blue
DR ROUND BUR	Ø 4.0	130	120.100.220	FMS blue
DR ROUND BUR	Ø 5.5	130	120.200.220	FMS blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver


Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	020.101.030	FMS
DR WAVE CUTTER	Ø 5.0	130	020.201.030	FMS
DR PIT BULL	Ø 3.0	130	020.301.040	FMS
DR PIT BULL	Ø 4.0	130	020.101.040	FMS
DR PIT BULL	Ø 5.0	130	020.201.040	FMS
DR PIT BULL	Ø 7.0	130	020.701.040	FMS
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	020.101.036	FMS
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	020.201.043	FMS
DR PIT BULL	Ø 4.0	155	020.101.037	FMS
DR PIT BULL	Ø 5.0	155	020.201.037	FMS
DR BENDING PIT BULL	Ø 3.0	130	020.301.140	FMS
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	020.101.041	FMS
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	020.101.042	FMS
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	020.101.050	FMS
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	020.101.051	FMS
DR BARREL BUR	Ø 4.0	130	130.100.020	FMS
DR BARREL BUR	Ø 5.5	130	130.200.020	FMS
DR ROUND BUR	Ø 4.0	130	120.100.020	FMS
DR ROUND BUR	Ø 5.5	130	120.200.020	FMS

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	020.101.528	FMS
DR WAVE CUTTER	Ø 5.0	130	020.201.528	FMS
DR PIT BULL	Ø 4.0	130	020.101.538	FMS
DR PIT BULL	Ø 5.0	130	020.201.538	FMS
DR PIT BULL	Ø 7.0	130	020.701.538	FMS
DR PIT BULL	Ø 4.0	155	020.101.535	FMS
DR PIT BULL	Ø 5.0	155	020.201.535	FMS
DR SYNOVIA PIT BULL	Ø 4.0	130	020.101.534	FMS
DR BARREL BUR	Ø 4.0	130	130.100.519	FMS
DR BARREL BUR	Ø 5.5	130	130.200.519	FMS
DR ROUND BUR	Ø 4.0	130	120.100.519	FMS
DR ROUND BUR	Ø 5.5	130	120.200.519	FMS

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES

**DYONICS
S&N**

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	021.101.230	DRSN blue
DR WAVE CUTTER	Ø 5.0	130	021.201.230	DRSN blue
DR PIT BULL	Ø 3.0	130	021.301.240	DRSN blue
DR PIT BULL	Ø 4.0	130	021.101.240	DRSN blue
DR PIT BULL	Ø 5.0	130	021.201.240	DRSN blue
DR PIT BULL	Ø 7.0	130	021.701.240	DRSN blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.236	DRSN blue
DR BENDING PIT BULL	Ø 3.0	130	021.301.340	DRSN blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.241	DRSN blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.242	DRSN blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.250	DRSN blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.251	DRSN blue
DR BARREL BUR	Ø 4.0	130	131.100.220	DRSN blue
DR BARREL BUR	Ø 5.5	130	131.200.220	DRSN blue
DR ROUND BUR	Ø 4.0	130	121.100.220	DRSN blue
DR ROUND BUR	Ø 5.5	130	121.200.220	DRSN blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	021.101.030	DRSN
DR WAVE CUTTER	Ø 5.0	130	021.201.030	DRSN
DR PIT BULL	Ø 3.0	130	021.301.040	DRSN
DR PIT BULL	Ø 4.0	130	021.101.040	DRSN
DR PIT BULL	Ø 5.0	130	021.201.040	DRSN
DR PIT BULL	Ø 7.0	130	021.701.040	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.036	DRSN
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	021.201.043	DRSN
DR PIT BULL	Ø 4.0	155	021.101.037	DRSN
DR PIT BULL	Ø 5.0	155	021.201.037	DRSN
DR BENDING PIT BULL	Ø 3.0	130	021.301.140	DRSN
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.041	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.042	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.050	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.051	DRSN
DR BARREL BUR	Ø 4.0	130	131.100.020	DRSN
DR BARREL BUR	Ø 5.5	130	131.200.020	DRSN
DR ROUND BUR	Ø 4.0	130	121.100.020	DRSN
DR ROUND BUR	Ø 5.5	130	121.200.020	DRSN

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	021.101.528	DRSN
DR WAVE CUTTER	Ø 5.0	130	021.201.528	DRSN
DR PIT BULL	Ø 4.0	130	021.101.538	DRSN
DR PIT BULL	Ø 5.0	130	021.201.538	DRSN
DR PIT BULL	Ø 7.0	130	021.701.538	DRSN
DR PIT BULL	Ø 4.0	155	021.101.535	DRSN
DR PIT BULL	Ø 5.0	155	021.201.535	DRSN
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.534	DRSN
DR BARREL BUR	Ø 4.0	130	131.100.519	DRSN
DR BARREL BUR	Ø 5.5	130	131.200.519	DRSN
DR ROUND BUR	Ø 4.0	130	121.100.519	DRSN
DR ROUND BUR	Ø 5.5	130	121.200.519	DRSN

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES **STRYKER**

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	022.101.230	DRST blue
DR WAVE CUTTER	Ø 5.0	130	022.201.230	DRST blue
DR PIT BULL	Ø 3.0	130	022.301.240	DRST blue
DR PIT BULL	Ø 4.0	130	022.101.240	DRST blue
DR PIT BULL	Ø 5.0	130	022.201.240	DRST blue
DR PIT BULL	Ø 7.0	130	022.701.240	DRST blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.236	DRST blue
DR BENDING PIT BULL	Ø 3.0	130	022.301.340	DRST blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.241	DRST blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.242	DRST blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.250	DRST blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.251	DRST blue
DR BARREL BUR	Ø 4.0	130	132.100.220	DRST blue
DR BARREL BUR	Ø 5.5	130	132.200.220	DRST blue
DR ROUND BUR	Ø 4.0	130	122.100.220	DRST blue
DR ROUND BUR	Ø 5.5	130	122.200.220	DRST blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	022.101.030	DRST
DR WAVE CUTTER	Ø 5.0	130	022.201.030	DRST
DR PIT BULL	Ø 3.0	130	022.301.040	DRST
DR PIT BULL	Ø 4.0	130	022.101.040	DRST
DR PIT BULL	Ø 5.0	130	022.201.040	DRST
DR PIT BULL	Ø 7.0	130	022.701.040	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.036	DRST
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	022.201.043	DRST
DR PIT BULL	Ø 4.0	155	022.101.037	DRST
DR PIT BULL	Ø 5.0	155	022.201.037	DRST
DR BENDING PIT BULL	Ø 3.0	130	022.301.140	DRST
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.041	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.042	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.050	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.051	DRST
DR BARREL BUR	Ø 4.0	130	132.100.020	DRST
DR BARREL BUR	Ø 5.5	130	132.200.020	DRST
DR ROUND BUR	Ø 4.0	130	122.100.020	DRST
DR ROUND BUR	Ø 5.5	130	122.200.020	DRST

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	022.101.528	DRST
DR WAVE CUTTER	Ø 5.0	130	022.201.528	DRST
DR PIT BULL	Ø 4.0	130	022.101.538	DRST
DR PIT BULL	Ø 5.0	130	022.201.538	DRST
DR PIT BULL	Ø 7.0	130	022.701.538	DRST
DR PIT BULL	Ø 4.0	155	022.101.535	DRST
DR PIT BULL	Ø 5.0	155	022.201.535	DRST
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.534	DRST
DR BARREL BUR	Ø 4.0	130	132.100.519	DRST
DR BARREL BUR	Ø 5.5	130	132.200.519	DRST
DR ROUND BUR	Ø 4.0	130	122.100.519	DRST
DR ROUND BUR	Ø 5.5	130	122.200.519	DRST

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES **ARTHREX**

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	023.101.230	DRAR blue
DR WAVE CUTTER	Ø 5.0	130	023.201.230	DRAR blue
DR PIT BULL	Ø 3.0	130	023.301.240	DRAR blue
DR PIT BULL	Ø 4.0	130	023.101.240	DRAR blue
DR PIT BULL	Ø 5.0	130	023.201.240	DRAR blue
DR PIT BULL	Ø 7.0	130	023.701.240	DRAR blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.236	DRAR blue
DR BENDING PIT BULL	Ø 3.0	130	023.301.340	DRAR blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.241	DRAR blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.242	DRAR blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.250	DRAR blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.251	DRAR blue
DR BARREL BUR	Ø 4.0	130	133.100.220	DRAR blue
DR BARREL BUR	Ø 5.5	130	133.200.220	DRAR blue
DR ROUND BUR	Ø 4.0	130	123.100.220	DRAR blue
DR ROUND BUR	Ø 5.5	130	123.200.220	DRAR blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	023.101.030	DRAR
DR WAVE CUTTER	Ø 5.0	130	023.201.030	DRAR
DR PIT BULL	Ø 3.0	130	023.301.040	DRAR
DR PIT BULL	Ø 4.0	130	023.101.040	DRAR
DR PIT BULL	Ø 5.0	130	023.201.040	DRAR
DR PIT BULL	Ø 7.0	130	023.701.040	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.036	DRAR
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	023.201.043	DRAR
DR PIT BULL	Ø 4.0	155	023.101.037	DRAR
DR PIT BULL	Ø 5.0	155	023.201.037	DRAR
DR BENDING PIT BULL	Ø 3.0	130	023.301.140	DRAR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.041	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.042	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.050	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.051	DRAR
DR BARREL BUR	Ø 4.0	130	133.100.020	DRAR
DR BARREL BUR	Ø 5.5	130	133.200.020	DRAR
DR ROUND BUR	Ø 4.0	130	123.100.020	DRAR
DR ROUND BUR	Ø 5.5	130	123.200.020	DRAR

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	023.101.528	DRAR
DR WAVE CUTTER	Ø 5.0	130	023.201.528	DRAR
DR PIT BULL	Ø 4.0	130	023.101.538	DRAR
DR PIT BULL	Ø 5.0	130	023.201.538	DRAR
DR PIT BULL	Ø 7.0	130	023.701.538	DRAR
DR PIT BULL	Ø 4.0	155	023.101.535	DRAR
DR PIT BULL	Ø 5.0	155	023.201.535	DRAR
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.534	DRAR
DR BARREL BUR	Ø 4.0	130	133.100.519	DRAR
DR BARREL BUR	Ø 5.5	130	133.200.519	DRAR
DR ROUND BUR	Ø 4.0	130	123.100.519	DRAR
DR ROUND BUR	Ø 5.5	130	123.200.519	DRAR

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES

WOLF

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	025.101.230	DRWO blue
DR WAVE CUTTER	Ø 5.0	130	025.201.230	DRWO blue
DR PIT BULL	Ø 3.0	130	025.301.240	DRWO blue
DR PIT BULL	Ø 4.0	130	025.101.240	DRWO blue
DR PIT BULL	Ø 5.0	130	025.201.240	DRWO blue
DR PIT BULL	Ø 7.0	130	025.701.240	DRWO blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.236	DRWO blue
DR BENDING PIT BULL	Ø 3.0	130	025.301.340	DRWO blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.241	DRWO blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.242	DRWO blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.250	DRWO blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.251	DRWO blue
DR BARREL BUR	Ø 4.0	130	135.100.220	DRWO blue
DR BARREL BUR	Ø 5.5	130	135.200.220	DRWO blue
DR ROUND BUR	Ø 4.0	130	125.100.220	DRWO blue
DR ROUND BUR	Ø 5.5	130	125.200.220	DRWO blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	025.101.030	DRWO
DR WAVE CUTTER	Ø 5.0	130	025.201.030	DRWO
DR PIT BULL	Ø 3.0	130	025.301.040	DRWO
DR PIT BULL	Ø 4.0	130	025.101.040	DRWO
DR PIT BULL	Ø 5.0	130	025.201.040	DRWO
DR PIT BULL	Ø 7.0	130	025.701.040	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.036	DRWO
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	025.201.043	DRWO
DR PIT BULL	Ø 4.0	155	025.101.037	DRWO
DR PIT BULL	Ø 5.0	155	025.201.037	DRWO
DR BENDING PIT BULL	Ø 3.0	130	025.301.140	DRWO
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.041	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.042	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.050	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.051	DRWO
DR BARREL BUR	Ø 4.0	130	135.100.020	DRWO
DR BARREL BUR	Ø 5.5	130	135.200.020	DRWO
DR ROUND BUR	Ø 4.0	130	125.100.020	DRWO
DR ROUND BUR	Ø 5.5	130	125.200.020	DRWO

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	025.101.528	DRWO
DR WAVE CUTTER	Ø 5.0	130	025.201.528	DRWO
DR PIT BULL	Ø 4.0	130	025.101.538	DRWO
DR PIT BULL	Ø 5.0	130	025.201.538	DRWO
DR PIT BULL	Ø 7.0	130	025.701.538	DRWO
DR PIT BULL	Ø 4.0	155	025.101.535	DRWO
DR PIT BULL	Ø 5.0	155	025.201.535	DRWO
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.534	DRWO
DR BARREL BUR	Ø 4.0	130	135.100.519	DRWO
DR BARREL BUR	Ø 5.5	130	135.200.519	DRWO
DR ROUND BUR	Ø 4.0	130	125.100.519	DRWO
DR ROUND BUR	Ø 5.5	130	125.200.519	DRWO

...the sharpest one on the market!



SHAVER BLADES & BURS FOR ARTHROSCOPY HANDPIECES

STORZ

Single use instruments into sterile packaging with 5 pieces inside

Plastic blue hub

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	026.101.230	DRSZ blue
DR WAVE CUTTER	Ø 5.0	130	026.201.230	DRSZ blue
DR PIT BULL	Ø 3.0	130	026.301.240	DRSZ blue
DR PIT BULL	Ø 4.0	130	026.101.240	DRSZ blue
DR PIT BULL	Ø 5.0	130	026.201.240	DRSZ blue
DR PIT BULL	Ø 7.0	130	026.701.240	DRSZ blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.236	DRSZ blue
DR BENDING PIT BULL	Ø 3.0	130	026.301.340	DRSZ blue
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.241	DRSZ blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.242	DRSZ blue
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.250	DRSZ blue
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.251	DRSZ blue
DR BARREL BUR	Ø 4.0	130	136.100.220	DRSZ blue
DR BARREL BUR	Ø 5.5	130	136.200.220	DRSZ blue
DR ROUND BUR	Ø 4.0	130	126.100.220	DRSZ blue
DR ROUND BUR	Ø 5.5	130	126.200.220	DRSZ blue

Single use instruments into sterile packaging with 5 pieces inside

Metal hub silver

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	026.101.030	DRSZ
DR WAVE CUTTER	Ø 5.0	130	026.201.030	DRSZ
DR PIT BULL	Ø 3.0	130	026.301.040	DRSZ
DR PIT BULL	Ø 4.0	130	026.101.040	DRSZ
DR PIT BULL	Ø 5.0	130	026.201.040	DRSZ
DR PIT BULL	Ø 7.0	130	026.701.040	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.036	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	026.201.043	DRSZ
DR PIT BULL	Ø 4.0	155	026.101.037	DRSZ
DR PIT BULL	Ø 5.0	155	026.201.037	DRSZ
DR BENDING PIT BULL	Ø 3.0	130	026.301.140	DRSZ
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.041	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.042	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.050	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.051	DRSZ
DR BARREL BUR	Ø 4.0	130	136.100.020	DRSZ
DR BARREL BUR	Ø 5.5	130	136.200.020	DRSZ
DR ROUND BUR	Ø 4.0	130	126.100.020	DRSZ
DR ROUND BUR	Ø 5.5	130	126.200.020	DRSZ

Reusable instruments into sterile packaging with 1 piece inside

Metal hub with gold surface

Description	Ø (mm)	L (mm)	No. Ref.	Hub
DR WAVE CUTTER	Ø 4.0	130	026.101.528	DRSZ
DR WAVE CUTTER	Ø 5.0	130	026.201.528	DRSZ
DR PIT BULL	Ø 4.0	130	026.101.538	DRSZ
DR PIT BULL	Ø 5.0	130	026.201.538	DRSZ
DR PIT BULL	Ø 7.0	130	026.701.538	DRSZ
DR PIT BULL	Ø 4.0	155	026.101.535	DRSZ
DR PIT BULL	Ø 5.0	155	026.201.535	DRSZ
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.534	DRSZ
DR BARREL BUR	Ø 4.0	130	136.100.519	DRSZ
DR BARREL BUR	Ø 5.5	130	136.200.519	DRSZ
DR ROUND BUR	Ø 4.0	130	126.100.519	DRSZ
DR ROUND BUR	Ø 5.5	130	126.200.519	DRSZ

...the sharpest one on the market!



DR BENDING PIT BULL® pat.

- The DR BENDING PIT BULL® pat. Ø3mm will be supplied straight and is bendable in every direction with the DR BENDING TOOL

Arthroscopy applications:

hip, knee, shoulder | elbow, ankle | wrist

DR BENDING TOOL 020.101.103

- reusable
- supplied non sterile



Instruction

- 3 different radius for optimal results
- introduce the DR BENDING PIT BULL in the DR BENDING TOOL
- bend it in your preferred angle position

DR REPROCESSING – SET

for reusable DR SHAVER BLADES & BURS



DR REPROCESSING - SET for hub	No. Ref.
DROR	140.000.220
FMS	140.000.220
DRSN Dyonics	140.000.220
DRAR Arthrex	140.000.220
DRST Stryker	140.000.260
DRWO Wolf	140.000.261
DRSZ Storz	140.000.262

Content DR REPROCESSING - SET

QTY	Description			No. Ref.
1	DR CASE			140.000.250
2	DR BOX M			140.000.251
1	DR PREWASH-TOOL			140.000.212
1	DR WASH-UNIT:			
	- DROR			140.000.213
	- FMS			140.000.213
	- DRSN Dyonics			140.000.213
	- DRAR Arthrex			140.000.213
	- DRST Stryker			140.000.218
	- DRWO Wolf			140.000.219
	- DRSZ Storz			140.000.259
1	DR LOUPE			140.000.209
1	DR HOT AIR GUN			140.000.210
1	DR GRAISSE (Tube 30g)			100.803.103
1	DR TUBE SET / 50 pieces	DR SHAVER	Ø4mm	140.000.255
1	DR TUBE SET / 50 pieces	DR SHAVER	Ø5mm	140.000.256
1	DR TUBE SET / 50 pieces	DR BUR	Ø4mm	140.000.257
1	DR TUBE SET / 50 pieces	DR BUR	Ø5.5mm	140.000.258




DR DISTRACTOR_{pat.} for hip arthroscopy

DR DISTRACTOR

Precision mechanism allowing optimal adaptation to the axis of the DR SCREW

Stability and power due to the angular stability of the DR DISTRACTOR and screw

Simple and effective control of distraction



ADVANTAGES FOR THE PATIENTS

No lesion on the nervus pudendus
(no impotence and erection problems)

No sciatic nerve lesion (no foot drop)

No femoral nerve lesion

No perineal lesion due to compression
on a traction table post


ADVANTAGES FOR THE SURGEON

No time pressure when performing reconstructive
and interventional hip arthroscopy


Simple handling and installation

Important and comfortable working area


Optimal DV-View, hinge mechanism avoiding the
shadowing effect while using the image intensifier



DR DISTRACTOR INSTRUMENTS



DR Instruments: Type S+M



DR Instruments: Type L

DR SCREW mod.dep.



DR SCREW for DR DISTRACTOR

DR SCREW K5 Cortical screw L175mm / type S+M

DR SCREW K6 Cortical screw L175mm / type S+M

DR SCREW K6 Cortical screw L210mm / type L

up to 60 kg body weight
more than 60 kg body weight
more than 60 kg body weight



Self cutting



Large 3-edge surfaces for optimum insertion
and removal of screws

DR BIOPSY – NEEDLE (for DR DISTRACTOR)



Extremely resistant and steady
Supplied sterile
Single use product

DR CKC – SET mod.dep.



The set comprises 2 parts:

- 1 highly-resistant arthroscopy knife with a long shaft
- 1 nitinol guide wire

Can also be used as an outflow cannula at the beginning of hip arthroscopy

Supplied sterile

Single use product



Applications

Hip Arthroscopy:

- Capsulotomy
- Labral detachment

Shoulder Arthroscopy:

- Biceps tenotomy

DR DISTRACTOR AND DR INSTRUMENTS

Description	sterile	non sterile	No. Ref.
DR DISTRACTOR incl. DR INSTRUMENTS TYPE S+M		x	100.803.088
DR INSTRUMENTS TYPE S+M		x	100.803.086
DR INSTRUMENTS TYPE L		x	100.803.085

Single use


DR SCREW K5 L175mm	packaging with 4 pieces inside	x		100.803.032
DR SCREW K6 L175mm	packaging with 4 pieces inside	x		100.803.033
DR SCREW K6 L210mm TYPE L	packaging with 4 pieces inside	x		100.803.025
DR BIOPSY - NEEDLE Ø2.1mm (for DR DISTRACTOR)	packaging with 2 pieces inside	x		100.803.111
DR CKC - SET Ø4mm L255mm	packaging with 5 pieces inside	x		100.803.202

Spare parts

DR SCREW NUT M6		x		100.803.030
DR GUIDING TUBES TYPE M		x		100.803.048
DR GUIDING TUBES TYPE S		x		100.803.049
DR GUIDING TUBES TYPE L		x		100.803.073
DR ADAPTER FOR DRILL		x		100.803.046
DR AO ADAPTER FOR DRILL		x		100.803.112
DR DRILL BIT Ø4mm for DR SCREW K5		x		100.803.091
DR DRILL BIT Ø4.5mm for DR SCREW K6		x		100.803.090
DR DRILL SOCKET K5		x		100.803.051
DR DRILL SOCKET K6		x		100.803.065
DR HOLDER FOR DRILL GUIDE		x		100.803.062
DR NAIL M+S K5		x		100.803.050
DR NAIL M+S K6		x		100.803.064
DR GRAISSE (Tube 30g)		x		100.803.103
DR MULTIPLE KEY		x		100.803.047
DR ALLEN KEY		x		100.803.061

DR HAI – SET

Hip Arthroscopy Instrument Set 100.803.101



DR HAI - SET (without 70° optics / L302mm)

100.803.101

Content DR HAI – SET

QTY	Description	No. Ref.
1	DR HANDLE	100.803.319
1	DR COUNTERBRANCE CANULATED	100.803.317
1	DR HOOK	100.803.306
2	DR COUNTERBRANCE	100.803.307
1	DR HALFTUBE	100.803.316
2	DR TUBE CONNECTOR	100.803.324
1	DR ARTHROSCOPY SHAFT (Ø6mm / with 2 stopcocks)	100.803.320
3	DR WORKING TUBE 70° / Ø7mm / Ø6mm	100.803.321
3	DR WORKING TUBE 70° / Ø5.5mm / Ø4.5mm	100.803.322
1	DR TUBE END PIECE	100.803.323
1	DR INSTRUMENT CASE	100.803.131
additional		
3	DR WORKING TUBE 30° / Ø7mm / Ø6mm	100.803.340
3	DR WORKING TUBE 30° / Ø5.5mm / Ø4.5mm	100.803.341

DR ARTHRO – TOWER





- MONITOR
- RECORDER
- CAMERA SYSTEM
- LIGHT SOURCE
 - COLD LIGHT CABLE
- UNIVERSAL PERISTALTIC PUMP FOR ARTHROSCOPY
- SHAVER SYSTEM
- ARTHROSCOPE
- ARTHROSCOPY TROCAR
- OBTURATOR

DR RACK

DR RACK for EQUIPMENT 200.000.043



DR MONITOR HOLDER 200.000.079
 horizontally and vertically adjustable

DR MONITOR 21" 200.000.042

DR CAMERA SET - HD Quality



DR CONSOLE for CAMERA 200.000.054

DR CAMERA 200.000.039
 with Vario Adapter for picture zoom



DR ARTHROSCOPE 30° / Ø4mm / L175mm (DR MEDICAL) 200.000.052

DR ARTHROSCOPE 70° / Ø4mm / L302mm (DR MEDICAL) 200.000.053

DR ARTHROSCOPE 30° / Ø4mm / L302mm (Baholzer) 100.803.410

DR ARTHROSCOPE 70° / Ø4mm / L302mm (Baholzer) 100.803.420



DR TROCAR (ARTHROSCOPY SHAFT)

Ø6mm / blue 200.000.055

DR OBTURATOR Ø 4.2mm / blue 200.000.056



DR CASE L290mm 200.000.216
 for DR ARTHROSCOPE L175mm

DR CASE L450mm 200.000.329
 for DR ARTHROSCOPE L302mm



DR RECORDER 200.000.041



DR KEYBOARD DE 200.000.074


DR KEYBOARD EN 200.000.075

DR KEYBOARD FR 200.000.076

DR LIGHT SOURCE



DR LED COLD LIGHT SOURCE 180W 200.000.031



DR CASE for DR HANDPIECE and DR LED COLD LIGHT CABLE 140.000.270
for Handpiece protection and washing process



DR LED COLD LIGHT CABLE 200.000.071


DR CONSOLE SET



DR CONSOLE for DR SHAVER 200.000.003



DR HANDPIECE 200.000.010
for DR SHAVER with DROR-hub
Milling- and Oscillation-mode up to 5000 rpm
Washing-Process History-File / long durability / always functional



DR INTERFACE CABLE 200.000.040
to Fluid pump for
FMS Duo+, Zimmer Biomet and Medical Vision



DR FOOTSWITCH
200.000.020

DR FLUID PUMP SET




DR UNIVERSAL FLUID PUMP
 200.000.034



DR REMOTE CONTROLLER
 200.000.035



DR DAY USE PUMP PART
 200.000.036



DR PATIENT PART SINGLE USE
 200.000.037

DR ARTHRO - TOWER Overview Equipment

DR RACK

Description	No. Ref.
DR RACK for EQUIPMENT	200.000.043
DR MONITOR HOLDER horizontally and vertically adjustable	200.000.079
DR MONITOR 21"	200.000.042

DR CAMERA SET - HD Quality (Light Source, Camera and Recorder)

Description	No. Ref.
DR CONSOLE for CAMERA	200.000.054
DR CAMERA WITH VARIO ADAPTER	200.000.039
DR ARTHROSCOPE 30° / Ø4mm / L175mm (DR MEDICAL)	200.000.052
DR ARTHROSCOPE 70° / Ø4mm / L302mm (DR MEDICAL)	200.000.053
DR ARTHROSCOPE 30° / Ø4mm / L302mm (Baholzer)	100.803.410
DR ARTHROSCOPE 70° / Ø4mm / L302mm (Baholzer)	100.803.420
DR TROCAR Ø6mm / blue	200.000.055
DR OBTURATOR Ø4.2mm / blue	200.000.056
DR CASE L290mm for DR ARTHROSCOPE L175mm	200.000.216
DR CASE L450mm for DR ARTHROSCOPE L302mm	200.000.329
DR RECORDER	200.000.041
DR KEYBOARD / DE	200.000.074
DR KEYBOARD / EN	200.000.075
DR KEYBOARD / FR	200.000.076

DR LIGHT SOURCE

Description	No. Ref.
DR LED COLD LIGHT SOURCE 180W	200.000.031
DR LED COLD LIGHT CABLE L3m / blue	200.000.071
DR CASE for DR HANDPIECE and DR LED COLD LIGHT CABLE	140.000.270

DR CONSOLE SET (Shaver Console and Handpiece for Arthroscopic intervention)

Description	No. Ref.
DR CONSOLE for DR SHAVER	200.000.003
DR HANDPIECE manual for DROR-hub	200.000.010
DR INTERFACE CABLE L3m (to Fluid pump for FMS Duo+, Zimmer Biomet and Medical Vision)	200.000.040
DR FOOTSWITCH with cable	200.000.020

DR FLUID PUMP SET (Arthroscopic Pump)

Description	No. Ref.
DR UNIVERSAL FLUID PUMP	200.000.034
DR REMOTE CONTROLLER	200.000.035
DR DAY USE PUMP PART	200.000.036
DR PATIENT PART SINGLE USE	200.000.037

design, development and production of equipment for health industry and hospitals
HIGH PRECISION INSTRUMENTS | SWISS MADE

SWITZERLAND
DR MEDICAL 
for your arthroscopy solution

DR MEDICAL AG
Ziegelmatteareal, Bau 6
CH-4503 Solothurn
SWITZERLAND

T. + 41 32 621 60 01
F. + 41 32 621 60 02
info@dr-medical.ch
www.dr-medical.ch
Version02_EN