



Non-slip “Step-on mat” for a safe footing of the medical staff



Non-slip backside

XL



Characteristics and benefits

- Safe footing for medical staff
- XL-format: 116 x 76 cm
- Individually cuttable
- Extra thin
- Non-slip backside
- Absorbs up to 1.7l of fluid

- + All surgery and examination areas in hospitals e.g. care and intensive care units, outpatient units, burn units
- + Nursing homes
- + Rest homes
- + Hospices
- + Doctors' practices
- + Home care services

Product information

Size	Absorption area	Packaging unit	Item No.
116 x 76 cm	116 x 76 cm	50 pcs	LMB11600760-622

SuperCore®: Patent held by the parent company of curea medical GmbH, McAiraid's Vliesstoffe GmbH

Picture credits: © curea medical GmbH

Disclaimer:

In spite of diligence the publisher cannot assume any liability for accuracy. We cannot assume any liability for information about indication and areas of application. Subject to modifications.

Protected brand names (trademarks) are not specially indicated. The absence of such notes may not be interpreted to mean that a trade name is unregistered. This leaflet and all of its content are protected by law. All rights reserved, especially the reprints (also of illustrations), any reproductions, lecture, radio, recordings and television broadcasts as well as storage in data-processing systems, even extracts or translations. Any use beyond the strict limits of the copyright law without the permission of curea medical GmbH is strictly prohibited and punishable by law.

© 2020 curea medical GmbH



curea medical
Innovation in Care

curea medical GmbH
Rudolf-Diesel-Straße 6
DE 37308 Heilbad Heiligenstadt
(Germany)

Telephone: +49 36071 9009500
Telefax: +49 36071 9009599
<http://www.curea-medical.de>
info@curea-medical.de

Registered Office:
Münsterstr. 61-65
DE 48565 Steinfurt (Germany)
AG Steinfurt, HRB 8735